









Setup



Sample Setup:

- Vial
- Aldehyde
- Stir bar
- No solvent

Machine Solvents

-  S1: Dichloromethane, anhydrous
-  S2: Hexafluoroisopropanol, anhydrous
-  S3: MeOH
-  S4: Diisopropylamine (175 mL) in THF (325 mL)
-  S5: –

Amount Aldehyde: 0-0.5 mmol

Solvent: Neat

Requirements: aldehyde soluble in CH₂Cl₂

Variants: Full sequence (11 hours) for highest conversion

Short sequence (5 hours) to obtain sample or
for electron poor aldehydes

Reaction Parameter Editing

Editable Parameters:


- Parameter 1: Imine formation temperature of cartridge (°C)
- Parameter 2: Imine formation temperature of reaction vial (°C)
- Parameter 3: Imine formation reaction time (seconds)
- Parameter 4: Cyclization temperature of cartridge (°C)
- Parameter 5: Cyclization temperature of reaction vial (°C)
- Parameter 6: Cyclization reaction time (seconds)

Editable Parts:


- Part 1: Substrate solvation step
- Part 2: Purification step

For more information see Application Note – N-Heterocycle

Additional Notes: 1) Check solvent level before setup. 2) aldehydes containing beta-heteroatom are not supported (catalyst complexation). 3) Run “Wash DCM” Sequence between each run.



Setup




Sample Setup:

- Vial
- Aldehyde / Ketone
- Amine
- Stir bar
- 5 mL solvent (see variants)


Machine Solvents

- S1: Dichloromethane, anhydrous
- S2: –
- S3: MeOH
- S4: Diisopropylamine (175 mL) in MeOH (325 mL)
- S5: –


Variants



Aldehyde + 1° Amine



Aldehyde + 2° Amine



Ketone + 1° Amine


Substrate ratios for best result:

Aldehyde : 1° Amine	Aldehyde : 2° Amine	Ketone : 1° Amine
1:2	2:1	1:1 (to 2:1)

Amount: 0-0.5 mmol

Solvent: CH₂Cl₂ : HFIP = 4:1 (5 mL) for best result

CH₂Cl₂ can be partially replaced with MeOH for solubility



Ketone + 2° Amine

Substrate ratio for best result: Ketone : 2° Amine = 2:1

Amount: 0-0.5 mmol

Solvent: Toluene (5 mL) – can be replaced with MeOH for solubility

Reaction Parameter Editing

Editable Parameters:

Ketones + 2° Amines:

- Parameter 1: Reduction step temperature cartridge
- Parameter 2: Reduction step temperature vial
- Parameter 3: Reaction time for reduction (sec.)
- Parameter 4: Amount of solvent for elution of product from SCX

Other sequences:

- Parameter 1: Reaction time for reduction (seconds)
- Parameter 2: Amount of solvent for elution of product from SCX

Editable Parts:

Aldehydes + 1° Amines

- Part 1: Excess amine scavenging
- Part 2: Purification step


Other sequences:

- Part 1: Purification step

Additional Notes: 1) Check solvent level before setup. 2) Run Wash DCM Sequence between each run.



Setup



Sample Setup:

- Vial
- Aldehyde / Ketone
- Stir bar
- 4 mL solvent (see right)

Amount Aldehyde/Ketone:	0-0.1 mmol
Solvent:	CH ₂ Cl ₂ : HFIP = 3:1 (4 mL)
Requirements:	aldehyde/ketone soluble in solvent mixture
Reaction time:	13 hours
For Cartridges:	P001-P025 and P101-P125 P001L-P025L and P101L-P125L

Reaction Parameter Editing

Editable Parameters:

- Parameter 1: Reaction time of reduction (seconds)
- Parameter 2: Amount of solvent for elution of product from SCX

For more information see Application Note – PROTAC (amine)

Editable Parts:

- Part 1: Purification step

Machine Solvents

- S1: Dichloromethane, anhydrous
- S2: –
- S3: MeOH
- S4: Diisopropylamine (175 mL) in Isopropanol (325 mL)
- S5: –

Additional Notes: 1) Avoid handling the CRBN products in protic solvents, due to instability of the Pomalidomide core ($t_{1/2}$ ~30 min in MeOH, 40°C). 2) Check solvent level before setup. 3) Run Wash DCM Sequence between each run.


Reaction sample preparation – PROTAC (active ester)

Version: 10.01.2022



Discontinued – use PROTAC via amide formation as replacement (P060-P085)

Setup



Sample Setup:

- Vial
- Primary or secondary amine
- Stir bar
- 4 mL solvent (see right)

Machine Solvents

- S1: Dichloromethane, anhydrous
- S2: –
- S3: Isopropanol
- S4: –
- S5: –

Amount Aldehyde/Ketone: 0-0.1 mmol

Solvent: 4 mL of preferably ACN, DMF, EtOAc, CH_2Cl_2

Requirements: amine soluble in chosen solvent

Variants: currently only free amine; amine salts possible soon

For Cartridges: P040-P065 and P140-P165

Reaction Parameter Editing

Editable Parameters:

- Parameter 1: Reaction time of amide formation

For more information see Application Note – PROTAC (amide)


Editable Parts:

- Part 1: Purification step

Additional Notes: 1) Avoid handling the CRBN products in protic solvents, due to instability of the Pomalidomide core ($t_{1/2}$ ~30 min in MeOH, 40°C). 2) Check solvent level before setup. 3) Run Wash DCM Sequence between each run.



Setup



Sample Setup:

- Vial
- Primary or secondary amine
- Stir bar
- 4 mL solvent (see right)

Machine Solvents

- S1: Dichloromethane, anhydrous
- S2: –
- S3: Isopropanol
- S4: –
- S5: –

Amount Aldehyde/Ketone: 0-0.1 mmol

Solvent: 2 mL DCM + 2mL MeCN

Requirements: amine soluble in chosen solvent

Variants: currently only free amine; amine salts possible soon

For Cartridges: P060-P085 and P160-P185

Reaction Parameter Editing

Editable Parameters:

- Parameter 1: Reaction time of amide formation

For more information see Application Note – PROTAC (amide)


Editable Parts:

- Part 1: Purification step

Additional Notes: 1) Avoid handling the CRBN products in protic solvents to long, due to instability of the Pomalidomide core ($t_{1/2}$ ~30 min in MeOH, 40°C). 2) Check solvent level before setup. 3) Run Wash DCM Sequence between each run.



Setup



Sample Setup:

- Vial
- Aldehyde / Ketone
- Stir bar
- 4 mL solvent (see right)

Amount Aldehyde/Ketone: 0-0.1 mmol

Solvent: HFIP (4 mL) for primary biotinylation reagents
TFE (4 mL) for secondary biotinylation reagents

Requirements: aldehyde/ketone soluble in solvent mixture

Reaction time: 15 hours

Reaction Parameter Editing

Editable Parameters:

- Parameter 1: Reaction time for reduction (seconds)

For more information see Application Note – Biotin (amine)

Editable Parts:

- Part 1: Purification step


Machine Solvents

- S1: Dichloromethane, anhydrous
- S2: –
- S3: MeOH
- S4: Diisopropylamine (175 mL) in MeOH (325 mL)
- S5: –

Additional Notes: 1) Check solvent level before setup. 2) Run Wash MeOH Sequence between each run to remove insoluble Biotin residues.



Setup



Sample Setup:

- Vial
- Amine or amine salt
- Stir bar
- 2.5 mL solvent (see right)

Amount Aldehyde/Ketone: 0-0.1 mmol

Solvent: DMF (2.5 mL)

Requirements: amine soluble in DMF

Reaction time: 6 hours

Reaction Parameter Editing

Editable Parameters:

- Parameter 1: Reaction time for amide formation (seconds)

For more information see Application Note – Biotin (amide)


Editable Parts:

- Part 1: Amine Freebasing
- Part 2: Purification step


Machine Solvents

- S1: Dichloromethane, anhydrous
- S2: –
- S3: MeOH
- S4: –
- S5: –

Additional Notes: 1) Check solvent level before setup. 2) Run Wash MeOH Sequence between each run to remove insoluble Biotin residues.



Setup



Sample Setup:

- Vial
- Amine
- Stir bar
- No solvent

Amount Amine: 0-0.5 mmol for B001
0-1.2 mmol for B002

Solvent: Neat

Requirements: amine MeOH soluble

Sequence: Choose between “free amine” and “amine salt” depending on the state of the amine substrate

Reaction Parameter Editing

Editable Parameters:

- Parameter 1: Reaction time for Boc protection (seconds)

For more information see Application Note – Boc protection


Editable Parts:

–

Machine Solvents

 S1: Dichloromethane, anhydrous

 S2: –

 S3: MeOH


 S4: –

 S5: –

Additional Notes: 1) Check solvent level before setup. 2) Run Wash DCM or MeOH Sequence between each run.



Setup



Sample Setup:

- Vial
- N-Boc starting material
- Stir bar
- No solvent

Amount Amine: 0-0.5 mmol

Solvent: Neat

Requirements: substrates soluble in Dimethoxyethane (DME)

Reaction Parameter Editing

Editable Parameters:

- Parameter 1: Reaction time for Boc protection (seconds)

For more information see Application Note – Boc deprotection


Editable Parts:

–

Machine Solvents

 S1: Dichloromethane, anhydrous


 S2: Dimethoxyethane (DME)

 S3: MeOH


 S4: –

 S5: –

Additional Notes: 1) Check solvent level before setup. 2) Run Wash DCM or MeOH Sequence between each run.



Setup



Sample Setup:

- Vial
- Alcohol and Nucleophile
- Stir bar
- No solvent

Machine Solvents

- S1: Dichloromethane, anhydrous
- S2: –
- S3: Ethanol
- S4: Diisopropylamine (175 mL) in Ethanol (325 mL)
- S5: –

Amount Substrates:

0.2 - 0.5 mmol

Substrate ratio:

1.0 equiv. alcohol / 1.5 equiv. nucleophile

Solvent:

Neat

(THF can be added for solubility)

Requirements:

starting materials soluble in CH₂Cl₂

Reaction time:

6.5 hours

Sequence:

standard Mitsunobu

Mitsunobu (basic substrates)

Current limitations and incompatibilities:

- Acid sensitive substrates not supported
- Substrates not soluble in CH₂Cl₂
- Yield may decrease up to 25% if less than 0.5 mmol of s.m. are used
- Not acidic nucleophiles (i.e. amides) cannot be removed by the purification system

Reaction Parameter Editing

Editable Parameters:


- Parameter 1: Reaction time (seconds)
- Parameter 2: Amount of solvent for elution of product from SCX

Editable Parts:


- Part 1: Purification step

For more information see Application Note – Mitsunobu

Additional Notes: 1) Check solvent level before setup. 2) Run Wash DCM Sequence between each run. If MeOH Wash Sequence was run before the reaction, some MeOH Addition product may be observed



Setup








Sample Setup:

- Vial
- Alcohol starting material
- Crossed stir bar (if available)
- No solvent

Amount Alcohol: 0-0.2 mmol

Solvent: Neat

Machine Solvents

-  S1: Dichloromethane, anhydrous
-  S2: Toluene, anhydrous with molecular sieves (4A)
-  S3: MeOH
-  S4: –
-  S5: Acetonitrile, HPLC grade

Reaction Parameter Editing

Editable Parameters:


- Parameter 1: Reaction time for Fluorination (seconds)

For more information see Application Note – Fluorination


Editable Parts:

- Part 1: Purification

Additional Notes: 1) Check solvent level before setup. 2) Run Wash DCM Sequence between each run. If MeOH wash was used, a DCM wash is recommended before the fluorination reaction, since residual MeOH could interfere with the reaction.



Setup



Sample Setup:

- Vial
- Primary Amine or primary amine salt
- Stir bar
- 1 mL solvent (see right side)
- KHCO₃ (see right side)

Amount Alcohol: 0-0.5 mmol

Solvent:

For primary amines:	Acetonitrile : H ₂ O = 1:1 (1.0 mL)
For primary amine salts:	Acetonitrile : H ₂ O = 1:1 (1.0 mL) + equimolar amount of solid KHCO ₃ to HX

Requirements: Soluble in Acetonitrile : H₂O mixture

Reaction Parameter Editing

Editable Parameters:

- Parameter 1: Reaction time for Azide Formation (seconds)

For more information see Application Note – Azide Formation

Editable Parts:

- Part 1: Purification





Machine Solvents

- S1: –
- S2: –
- S3: MeOH
- S4: –
- S5: Acetonitrile : H₂O = 1 : 1


Additional Notes: 1) Check solvent level before setup. 2) Run MeOH Wash Sequence between each run. 3) When using materials containing basic functional groups (e.g. trialkylamines, imidazoles, pyridines, anilines) the purification step needs to be disabled to avoid the product being scavenged.

Reaction sample preparation – Amide formation

Version: 10.01.2022








Setup



Sample Setup:

- Vial
- Amine
- Carboxylic acid
- Stir bar
- 4.0 mL solvent (see right)

Machine Solvents

-  S1: Dichloromethane, anhydrous
-  S2: –
-  S3: MeOH
-  S4: –
-  S5: –

Amount Amine:
Ratio Amine : Acid

0-0.5 mmol
1 : 1.1

Solvent:

4.0 mL (2.0 mL anhydrous EtOH + 2.0 mL anhydrous CH₂Cl₂)

Requirements:

substrates soluble in EtOH + CH₂Cl₂ mixture

Reaction time:

6 hours for secondary amines. Can be shortened to 3 hours for primary amines (see application note for examples)

Reaction Parameter Editing

Editable Parameters:


- Parameter 1: Reaction time for amide formation (seconds)

For more information see Application Note – Amide formation


Editable Parts:

- Part 1: Purification step

Additional Notes: 1) Check solvent level before setup. 2) Run MeOH wash sequence after reaction.



Setup



Sample Setup:

- Vial
- Silyl protected compound
- Stir bar
- 2.0 mL solvent (see right)

Machine Solvents

- S1: Dichloromethane, anhydrous
- S2: –
- S3: MeOH
- S4: –
- S5: –

Amount Substrate:

0-0.5 mmol
[Si] = TBS, TES, TIPS

Solvent:

- 2.0 mL MeOH for best result
- alternatively 2.0 mL EtOH or *i*PrOH (requires longer reaction time, e.g. 13 hours).
- in case of poor solubility a mixture of 1.0 mL CH₂Cl₂ + 1.0 mL MeOH can be used

Requirements:

substrates soluble in the solvent mixture

Reaction Parameter Editing

Editable Parameters:

- Parameter 1: Reaction time for deprotection reaction (seconds)

For more information see Application Note – Silyl deprotection

Editable Parts:

- Part 1: Purification step

Additional Notes: 1) Check solvent level before setup. 2) Run MeOH wash sequence between each run.